

Contributed by: Umair Khan Academy www.umairkhanacademy.com

F.Sc. Practical Syllabus

1. Short experiments

- a. Crystallization of Benzoic acid from water.
- b. Purify commercial preparation of NaCl by passing HCl gas.
- c. Separate the mixture of various inks by paper chromatography
- d. Separation and identification of lead and cadmium ions from the mixture by paper chromatography.
- e. Determine the heat of neutralization of NaOH and HCl.

2. Titration (volumetric analysis)

- a. Acid base titration
- b. Redox titration
 - i. $KMnO_4$ titration
 - ii. Iodimetry

3. Detections

- a. Detection of elements C, H, N, S and halogens in organic compounds.
- b. Detection of functional group. (Carboxylic acid, phenol and aldehyde are in notes)

4. Preparations

- a. Preparation of aspirin (acetyl salicylic acid)
- b. Preparation of iodoform
- c. Preparation of phenyl glucasazone.
- d. Preparation of copper ammine complex.

5. Inorganic analysis (qualitative analysis or salt analysis)

- a. Detection of Acid radical
- b. Detection of Basic radical

6. Estimation of barium in the given sample solution, gravimetrically as chromate

i

INSTRUCTION NOTE

- 1. In first 20 mints you have to write two experiments. Titration is of much importance. There are many experiments in practical notebook. Few important examples are included in notes. You have to write, same as mentioned in notes with supposed observations and calculation in exams, even writing in you practical notebook.
- 2. 2nd most important experiment is salt analysis (qualitative analysis), they can be performed with the help of scheme included in the notes.
- 3. Short experiments can be consulted in initial experiments from notebook.
- 4. If you want to learn more about experiments, such as
 - a. How to perform basic experiment.
 - b. How to make solutions.
 - c. What are the basic concepts of titration
 - d. How to perform salt analysis and how to write in exam.
 - e. Basic definitions and concepts.
 - f. Online help Then consult UMAIR KHAN ACADEMY (Channel) At YouTube Or Facebook page
 https://www.facebook.com/umairkhanacademy/

Some links are given below

1. Prepare 1 molar, 1000 ml solution of NaOH.

https://youtu.be/_q7Q4md-ZCY

2. Prepare 1 molar, 100 ml solution of NaOH.

https://youtu.be/5jfpd-4w9tI

3. Prepare 0.05 molar, 1000 ml solution of NaOH.

https://youtu.be/-pX3jhG7tbg

4. Prepare 1 molar, 1000 ml solution of oxalic Acid.

https://youtu.be/oyXtSfhAmPU

5. **Prepare 1 molar, 50 ml solution of H2SO4 from stock molarity (18M)** https://youtu.be/vQPOtyjuZxk

6. **Prepare 1 molar, 1000 ml solution of HCl from stock molarity (12M)** https://youtu.be/xXOtRmtB4Wc

7. Prepare 1 molar, 1000 ml solution of HCl from stock solution which is 37% of HCl having density 1.18 g/cm3.

https://youtu.be/Cnx4gt3wUZ0

8. **A complete guide to titration (Very much informative)** https://youtu.be/cD_GMCQd7UE

Quick Points to Learn

N-4

Some important acid base reactions

2. $2HCl + Na_2CO_3 \rightarrow H_2CO_3 + 2NaCl$	molar ratio 2:1
3. $H_2SO_4 + 2NaOH \rightarrow Na_2SO_4 + 2H_2O$	Molar ratio 1:2
4. $H_2SO_4 + 2KOH \rightarrow K_2SO_4 + 2H_2O$	Molar ratio 1:2
5. $(COOH)_2 + 2NaOH \rightarrow (COONa)_2 + 2H_2O$	Molar ratio 1:2
6. $H_2SO_4 + Na_2CO_3 \rightarrow Na_2SO_4 + H_2CO_3$ 7. $H_2SO_4 + 2NaHCO_3 \rightarrow Na_2SO_4 + H_2CO_3$	Molar ratio 1:1 Molar ratio 1:2
7. $H_2SO_4 + 2NaHCO_3 \rightarrow Na_2SO_4 + H_2CO_3$ 8. $CH_3COOH + NaOH \rightarrow CH_3COONa + H_2O$	
5 5 2	
<u>Indicator</u>	

- **1.** Phenolphthalein is used in case of strong bases such as NaOH, KOH.
- 2. Methyl orange is used in case of weak bases such as Na₂CO₃, Soap(NaOH)
- 3. <u>KMnO₄ Titration (KMnO₄ is used as oxidizing agent with following reducing reagents)</u>
- i. FeSO₄
- ii. FeSO₄.7H₂O
- iii. $FeSO_4.(NH_4)_2SO_4.6H_2O$
- **iv.** (COOH)₂
- **v.** (COOH)₂.2H₂O
- vi. $(COONa)_2$
- **vii.** (COOK)₂
- viii. (COONH₄)₂

Molar Wt. 152 g/mol Molar Wt. 278 g/mol Molar Wt. 392 g/mol (Mohar's salt) Molar Wt. 90 g/mol (Oxalic Acid) Molar Wt. 126 g/mol Molar Wt. 134 g/mol (Sod. oxalate) Molar Wt. 166 g/mol (Pot. oxalate) Molar Wt. 124 g/mol (Amm. oxalate)

Important redox reactions

- 1. $2KMnO_4 + 8H_2SO_4 + 10FeSO_4 .7H_2O \rightarrow K_2SO_4 + 2MnSO_4 + 5Fe_2(SO_4)_3 + 78H_2O$
- 2. $2KMnO_4 + 10FeSO_4(NH_4)_2SO_4.6H_2O + 8H_2SO_4 \rightarrow K_2SO_4 + 2MnSO_4 + 5Fe_2(SO_4)_3 + 10(NH_4)_2SO_4 + 68H_2O_4 + 68$
- **3.** $2KMnO_4 + 3H_2SO_4 + 5(COOH)_2.2H_2O \rightarrow K_2SO_4 + 2MnSO_4 + 10CO_2 + 18H_2O_2$
- 4. $2KMnO_4 + 8H_2SO_4 + 5(COONa)_2 \rightarrow K_2SO_4 + 2MnSO_4 + 10CO_2 + 8H_2O + 5Na_2SO_4$
- 5. $2KMnO_4 + 8H_2SO_4 + 5(COONH_4)_2 \rightarrow K_2SO_4 + 2MnSO_4 + 10CO_2 + 8H_2O + 5(NH4)_2SO_4$
- 6. $2KMnO_4 + 8H_2SO_4 + 5(COOK)_2 \rightarrow 6K_2SO_4 + 2MnSO_4 + 10CO_2 + 8H_2O_3 + 8H_2$
- 7. $2Na_2S_2O_3 + I_2 \rightarrow Na_2S_2O_6 + 2NaI$ (Iodimetry, redox reaction)

Standard Solution

The solution with known concentration (molarity) is called standard solution. **Molar Ratio**

Molar ratio is indicated by balanced chemical equation.

SHORT EXPERIMENTS

F.Sc Practical

Page No. 1

Question No. 1: Separate and identify the mixture of various inks by paper chromatography.

Apparatus and Chemicals

Whatman filter paper strip, chromatographic tank, beaker, given mixture of inks, solvent (Ethanol, water), capillary tube jet, scale, lead pencil

Procedure

- 1. Take a strip of Whatman filter paper, 20cm long.
- 2. Draw a line 2.5cm apart from one end of this strip with the help of lead pencil and mark 'X' in middle.
- 3. Put a small drop of ink mixture on the mark.
- 4. Now hang down the strip in chromatographic tank or beaker to dip one end (having mixture spot) in the solvent in a way, mixture spot remains outside from solvent.
- 5. Allow to stand the strip, till solvent travels ³/₄ part of it and different ink colour spots are visible.
- 6. Now, remove the strip, stand to dry.
- 7. Note the distance covered by solvent and different inks.
- 8. Calculate Rf value for each component by the given formula.

$R_f = \frac{Distance \; (cm) \; travelled \; by \; spot \; from \; original \; ilne}{Distance \; (cm) \; travelled \; by \; solvent \; from \; original \; ilne}$

<u>Result:</u> I separated the colours of different inks from given inks mixture and shown to the examiner.

Question No. 2: Separate and identify the mixture of Pb²⁺ and Cd²⁺ ions from their solutions by paper chromatography.

Apparatus and Chemicals

Whatman filter paper strip, chromatographic tank, beaker, given mixture of Pb^{2+} and Cd^{2+} ions, solvent (Acetone 2volume, *ter*-Butyl alcohol 2 volume, dil HNO₃ 1 volume), H₂S solution, capillary tube jet, scale, lead pencil.

Procedure

- 1. Take a strip of Whatman filter paper, 20cm long.
- 2. Draw a line 2.5cm apart from one end of this strip with the help of lead pencil and mark 'X' in middle.
- 3. Put a small drop of Pb^{2+} and Cd^{2+} mixture on the mark.
- 4. Now hang down the strip in chromatographic tank or beaker to dip one end (having mixture spot) in the solvent in a way, mixture spot remains outside from solvent.
- 5. Allow to stand the strip, till solvent travels $\frac{3}{4}$ part of it.
- 6. Now, remove the strip, stand to dry.
- Use locating agent (H₂S solution or ammonium sulphide soulution) to locate the spots of Pb²⁺ and Cd²⁺ ions.
- 8. Note the distance covered by solvent and Pb^{2+} and Cd^{2+} ions.
- 9. Calculate Rf value for each ion by the given formula.

SHORT EXPERIMENTS

F.Sc Practical

Page No. 2

 $R_{f} = \frac{Distance (cm) travelled by spot from original ilne}{Distance (cm) travelled by solvent from original ilne}$

<u>Result:</u> I separated the Pb²⁺ and Cd²⁺ ions from given mixture and shown to the examiner.

Question No. 3: Purify the given sample of benzoic acid by crystallization method.

Apparatus and Chemicals

Benzoic acid, water, spirit lamp, filters paper, glass rod, beaker, funnel, tripod stand,

Procedure

- 1. Take 50 ml water in 250 ml beaker and boil the water on tripod stand.
- 2. Add given sample of benzoic acid in small fractions with constant stirring till appropriate amount is added. (Do not make saturated or super saturated solutions)
- 3. Filter the hot solution through filter paper and funnel and filtrate is obtained in china dish.
- 4. Cool the filtrate, pure colourless, odourless crystalline needles of benzoic acid are obtained.
- 5. Dry the crystals between two folds of filter paper.

<u>Result:</u> I purified the given sample of benzoic acid by crystallization and shown to examiner.

Question No. 4: Purify the given sample of common salt by passing HCl gas.

Apparatus and Chemicals

Beaker, funnel, round bottom flask, gas burner or spirit lamp, glass tubing, NaCl, H₂SO₄

Procedure

- 1. Dissolve some quantity of impure salt in a beaker. Filter it in another beaker with filtering apparatus.
- 2. In a round bottom flask NaCl is reacted with Conc. H₂SO₄, dropped from thistle funnel, as a result HCl gas is formed which is passed through rubber tubing joined to the inverted funnel onto impure solution (filtrate of point 1) of sodium chloride.

 $2NaCl + H_2SO_4 \longrightarrow Na_2SO_4 + 2HCl$

- 3. White, pure crystals of sodium chloride are separated out due to common ion effect.
- 4. Dry them in folding of filter paper.

<u>Result:</u> I purified the given sample of NaCl and shown to examiner.

SHORT EXPERIMENTS

F.Sc Practical

UMAIR KHAN ACADEMY www.youtube.com/umairkhanacademy www. UmairKhanAcademy.com

Question No. 5: Determine the heat of neutralization of strong acid (HCl) and strong Base (NaOH).

Apparatus and Chemicals

Calorimeter with stirrer, thermometer, balance, beaker, glass rod, 0.1M HCl and 0.1M NaOH

Procedure

- 1. Take 50ml 0.1M NaOH solution in calorimeter and note the initial temperature (T1) by thermometer.
- 2. Now, add 50ml 0.1M HCl solution in calorimeter, neutralization reaction get started and heat is evolved.
- 3. Note the increasing temperature after each 10 seconds interval.
- 4. When temperature becomes uniform, it is noted as final temperature (T2)
- 5. Note the weight of calorimeter before adding solutions (M1) and after adding both solutions (M2).
- 6. Calculate the heat of neutralization as follows.

Volume of HCl or NaOH	=	V	=	50ml	
Molarity of HCl and NaOH	=	Μ	=	0.1M	
Specific Heat of		01		0.200 1/ 1/	
Calorimeter		S1	=	0.380 J/gK	
Specific Heat of Solution	=	S2	=	4.18 J/gK	
Rise in temperature	-	θ	=	T2- T1	
**			=	35.72-35	
			=	0.72 °C	
Heat of neutralization			=	$\frac{(M1S1+M2S2)}{1000}\times\frac{1000}{V\times M}$	
			_	$\underbrace{(50.5 \times 0.38 + 91.5 \times 4.18)}_{\times} \underbrace{1000}_{\times}$	
		_	$\frac{1000}{50 \times 0.1}$		
			=	57.4 kJ/mole	

<u>Result:</u> I determined heat of neutralization of strong acid (HCl) and strong base (NaOH) by calorimeter and I found it as -57.4 kJ/mole. (exothermic)

Page No. 3

Practical Notes (Titration) Page No. 4 **F.Sc. 2nd year Chemistry**

Question: 01

Question. 01					
<u>The given solution contains 6g of a mixture of</u>					
<u>HCl and NaCl dissolve per litre. Determine the</u>					
percentage compositi	percentage composition of the mixture. You are				
provided with 0.1M N	NaOH.				
Principle	It is an acid base titration.				
Standard Solution	0.1M NaOH				
Chemical Equation	$NaOH + HCl \rightarrow NaCl +$				
H ₂ O					
<u>Molar ratio</u>	$n_1 = 1, n_2 = 1$				
Indicator	Phenolphthalein				
<u>End point</u>	Light pink				
Procedure					

The given solution of HCl mixture is taken in burette and run down drop wise in the conical flask, having 10ml NaOH solution with few drops of phenolphthalein, till the end point is reached which is light pink. Three concordant readings are taken for volume of acid used. Observations

Initial reading 0 10 20 0 ns Acid $\frac{M_1V_1}{n_1}$ $\frac{M_1 10}{1}$ M_1	Final rea 10 20 30 Concord = = =		Volume of Acid used 10 ml 10 ml 10 ml eading (V ₁) = 10 m Base $\frac{M_2V_2}{n_2}$ 0.1 X 10 1 0.1 X 10 X 1		Observat	Initial reading 0 10 20 ions Acid $\frac{M_1V_1}{n_1}$
$\begin{array}{r} 0 \\ 10 \\ 20 \end{array}$ $\begin{array}{r} \text{Ons} \\ \text{Acid} \\ \frac{M_1 V_1}{n_1} \\ \frac{M_1 10}{1} \end{array}$	20 30 Concord	dant R	$\frac{10 \text{ ml}}{10 \text{ ml}}$ $\frac{10 \text{ ml}}{10 \text{ ml}}$ Reading (V ₁) = 10 m $\frac{\text{Base}}{\frac{M_2 V_2}{n_2}}$ $\frac{0.1 X 10}{1}$		$\frac{1}{2}$	$\frac{reading}{0}$ 10 20 $\frac{1}{10}$ $\frac{1}{20}$ $\frac{1}{10}$ $\frac{1}{10}$ $\frac{1}{10}$ $\frac{1}{10}$
$ \begin{array}{r} 10 \\ 20 \\ \hline M_1V_1 \\ \hline M_1 10 \\ \hline 1 \end{array} $	20 30 Concord	dant R	$\frac{10 \text{ ml}}{10 \text{ ml}}$ Beading (V ₁) = 10 m $\frac{Base}{\frac{M_2V_2}{n_2}}$ $\frac{0.1 X 10}{1}$		$\frac{1}{2}$	$ \begin{array}{r} 0 \\ 10 \\ 20 \\ \hline \text{ions} \\ \text{Acid} \\ \frac{M_1 V_1}{n_1} \\ \end{array} $
20 $\frac{\text{Dns}}{\text{Acid}}$ $\frac{M_1V_1}{n_1}$ $\frac{M_1 \ 10}{1}$	30 Concord	dant R	$\frac{10 \text{ ml}}{\text{leading } (V_1) = 10 \text{ m}}$ $\frac{\text{Base}}{\frac{M_2 V_2}{n_2}}$ $\frac{0.1 X 10}{1}$		2 3	$ \begin{array}{r} 0 \\ 10 \\ 20 \\ \hline \text{ions} \\ \text{Acid} \\ \frac{M_1 V_1}{n_1} \\ \end{array} $
$\frac{DONS}{Acid}$ $\frac{M_1V_1}{n_1}$ $\frac{M_1 \ 10}{1}$	Concore	dant R	Reading (V ₁) = 10 m Base $\frac{M_2V_2}{n_2}$ $\frac{0.1 \times 10}{1}$		3	$\frac{20}{\text{Acid}}$ $\frac{M_1V_1}{n_1}$
Acid $\frac{M_1V_1}{n_1}$ $\frac{M_1 \ 10}{1}$		dant R	Base $\frac{M_2V_2}{n_2}$ $\frac{0.1 \times 10}{1}$			$\frac{\text{Acid}}{\frac{M_1V_1}{n_1}}$
Acid $\frac{M_1V_1}{n_1}$ $\frac{M_1 \ 10}{1}$	=		$\frac{\underline{M}_2 V_2}{\underline{n}_2}$ $\frac{\underline{0.1 \times 10}}{1}$)	<u>Calculat</u>	Acid $\frac{M_1V_1}{n_1}$
$\frac{\frac{M_1V_1}{n_1}}{\frac{M_1 10}{1}}$	=		$\frac{\underline{M}_2 V_2}{\underline{n}_2}$ $\frac{\underline{0.1 \times 10}}{1}$		<u>Calculati</u>	Acid $\frac{M_1V_1}{n_1}$
$ \overline{n_1} \overline{n_1 \ 10} \overline{1} $	=		$\frac{\overline{n_2}}{0.1 \times 10}$			$\frac{M_1V_1}{n_1}$
$\frac{M_1 10}{1}$	-		$\frac{0.1 \times 10}{1}$			n_1
	- 0	<	1			_
	Ō	\mathbf{X}	_			M 10
M_1	Ō		0 1 V 10 V 1			<i>M</i> ₁ 10
			U. I A IU A I			1
			1 X 10			\mathbf{M}_{1}
arity of HCl	=		0.1 M			
lissolve per litr	e =	Mola	rity x Molecular. W	/ t.	Molarit	y of oxalic acid
					Amount	dissolve per
\sim	=		0.1 x 36.5		litre of o	xalic acid
	=		3.65g			
lissolve per litr	e =		6-3.65			
					Amount	dissolve per
	=		2.35g			odium oxalate
mposition of	=		3.65			
			<u> </u>		% age co	omposition of
	=		60.83%		oxalic ac	
mposition of	=		2.35			
-			$\frac{1}{6}$ X 100		% age co	omposition of
	=		39.17%		sodium o	-
	lissolve per litr	= lissolve per litre = = mposition of = mposition of =	= = lissolve per litre = = nposition of = = nposition of =	$= 0.1 \times 36.5$ $= 3.65g$ lissolve per litre = 6-3.65 $= 2.35g$ mposition of = $\frac{3.65}{6} \times 100$ $= 60.83\%$ mposition of = $\frac{2.35}{6} \times 100$	$= 0.1 \times 36.5$ $= 3.65g$ lissolve per litre = 6-3.65 $= 2.35g$ mposition of = $\frac{3.65}{6} \times 100$ $= 60.83\%$ mposition of = $\frac{2.35}{6} \times 100$	Amount $= 0.1 \times 36.5$ $= 3.65g$ $= 3.65g$ $= 6-3.65$ Amount = 2.35g $= 2.35g$ $= 60.83%$ Amount = 2.35g $= 60.83%$ Amount = 60.83%

UMAIR KHAN ACADEMY www.youtube.com/umairkhanacademy

www.youtube.com/umairkhanacadem www. UmairKhanAcademy.com

Question: 02

The given solution c	<u>ontains 8g of a mixture of</u>
Oxalic acid and sodi	<u>um oxalate dissolve in</u>
<u>1000ml. Determine t</u>	he percentage composition of
sodium oxalate the n	<u>nixture.</u>
Principle	It is an acid base titration.
Standard Solution	0.1M NaOH
Chemical Equation	
(COOH) ₂ + 2NaOH	\rightarrow (COONa) ₂ + 2H ₂ O
<u>Molar ratio</u>	(Acid) $n_1 = 1$, $n_2 = 2$
Indicator	Phenolphthalein
End point	Light pink
Procedure	

The given solution of oxalic acid mixture is taken in burette and run down drop wise in the conical flask, having 10ml NaOH solution with few drops of phenolphthalein, till the end point is reached which is light pink. Three concordant readings are taken for volume of acid used.

Observations Sr. NO Initial **Final reading** Volume of Acid reading used 10 0 10 ml 1 2 10 20 10 ml 3 20 30 10 ml Concordant Reading $(V_1) = 10$ ml **Calculations** Acid Base M_1V_1 M_2V_2 = n_1 n_2 $M_1 \, 10$ 0.1X10 = 2 1 0.1X10X1 M_1 = 2 X 10 Molarity of oxalic acid 0.05 M = Amount dissolve per Molarity x Molecular. =

Wt.

=

=

=

=

=

=

=

=

0.05 x 126

6.3 g

8-6.3

1.7g

 $\frac{6.3}{8}$ X 100

78.75%

 $\frac{1.7}{8}$ X 100

21.25%

Page No. 5 **F.Sc. 2nd year Chemistry**

Question: 03

The given solution c	The given solution contains 6.3g of (COOH)2.				
<u>xH₂O dissolve in litre (1000ml). Find out the</u>					
value of "x". You are provided with 0.1M NaOH.					
Principle	It is an acid base titration.				
Standard Solution	0.1M NaOH				
Chemical Equation					
$\overline{(\text{COOH})_{2.}\text{xH}_{2}\text{O}+2\text{Na}}$	$(COOH)_{2.x}H_{2}O + 2NaOH \rightarrow (COONa)_{2} + (2+x)H_{2}O$				
<u>Molar ratio</u>	(Acid) $n_1 = 1$, $n_2 = 2$				
Indicator	Phenolphthalein				
End point	Light pink				
Procedure					

The given solution of oxalic acid is taken in burette and run down drop wise in the conical flask, having 10ml NaOH solution with few drops of phenolphthalein, till the end point is reached which is light pink. Three concordant readings are taken for volume of acid used. <u>Observations</u>

Sr. NO	Initial	Final reading	Volume of Acid
	reading		used
1	0	10	10 ml
2	10	20	10 ml
3	20	30	10 ml
		Concordant l	Reading (V ₁) = 10 ml
<u>Calculat</u>	<u>ions</u>		
	Acid		Base
	M_1V_1	=	M_2V_2
	<i>n</i> ₁		n ₂
	$\frac{M_1 10}{1}$	=	<u>0.1 X 10</u>
	1		2
	\mathbf{M}_{1}	=	0. 1 X 10 X 1
			2 X 10
Molarit	y of oxalic acid		0.05 M
Amount	dissolve per	= Molar	rity x Molecular.
litre of o	xalic acid	Wt.	-
	6.3g/l	= 0.	05 X $(90 + 18x)$
	6.3	=	90 + 18x
	0.05		
	126	=	90 + 18 x
	126 -90	=	18x
	36		18x
	36	=	X
	18		
	2	=	Х
_			

<u>Result</u>

Water molecules of crystallization in oxalic acid are two

UMAIR KHAN ACADEMY www.youtube.com/umairkhanacademy www. UmairKhanAcademy.com

Question: 04 Determine the solubility of o

Determine the solubility of oxalic acid at room					
temperature. You are provided with 0.1M NaOH.					
<u>Principle</u>	It is an acid base titration.				
Standard Solution	0.1M NaOH				
Chemical Equation					
(COOH) _{2.} + 2NaOH -	\rightarrow (COONa) ₂ + 2H ₂ O				
<u>Molar ratio</u>	(Acid) $n_1 = 1$, $n_2 = 2$				
Indicator	Phenolphthalein				
<u>End point</u>	Light pink 🔹				
Procedure					

Clear saturated solution of oxalic acid at room temperature is prepared and 5ml of this solution is taken to make it 100 ml by diluting. This solution of oxalic acid is taken in burette and run down drop wise in the conical flask, having 10ml NaOH solution with few drops of phenolphthalein, till the end point is reached which is light pink. Three concordant readings are taken for volume of acid used.

Observations

<u>Observa</u>	<u>tions</u>			
Sr. NO	Initial	Final re	ading	Volume of Acid
	reading			used
1	0	10		10 ml
2	10	20		10 ml
3	20	30		10 ml
		Conco	rdant]	Reading $(V_1) = 10 \text{ ml}$
Calculat	<u>ions</u>			
	Acid			Base
	M_1V_1	=		M_2V_2
	n_1			n_2
	<i>M</i> ₁ 10	=		0.1X10
	1			2
	M_1	=	-	0.1 <i>X</i> 10 <i>X</i> 1
				2 X 10
Molaı	rity of oxalic	=		0.05 M
	acid			
Amoun	t dissolve per	r =	Mola	arity x
litre of	oxalic acid		Mole	ecular. Wt.
		=		0.05 x 126
		=		6.3 g/l
1000	ml of dilute	=		6.3g
oxalic	acid contains	5		_
1ml of	dilute oxalic	=		6.3
acio	d contains			1000
100 r	nl of dilute	=		$\frac{6.3}{1000}$ X 100
oxalic	acid contains	5		1000

Practical Notes (Titration) F.Sc. 2nd year Chemistry

Page No. 6

	=	0.63g
5ml of oxalic acid saturated solution	=	100 ml of dilute solution
100 ml of dilute solution contains	=	0.63g
5ml of oxalic acid solution contains	=	0.63g
1ml of oxalic acid solution	=	$\frac{0.63}{5}$
100ml of oxalic acid solution contains	=	$\frac{0.63}{5}$ X 100
	=	12.6g
D14		

Result

Solubility of oxalic acid at room temperature is 12.6g

UMAIR KHAN ACADEMY www.youtube.com/umairkhanacademy

www. UmairKhanAcademy.com

Question: 05				
<u>Determine the weight of monovalent metal M.</u>				
<u>When 4g of its hydroxide (MOH) are dissolve per</u>				
<u>litre in the given solution. You are provided with</u>				
<u>0.1M HCl.</u>				
Principle	It is an acid base titration.			
Standard Solution	0.1M HCl			
Chemical Equation	$MOH + HCl \rightarrow MCl + H_2O$			
<u>Molar ratio</u>	$n_1 = 1, n_2 = 1$			
Indicator	Phenolphthalein			
End point	Light pink			
Procedure				

..

^-

The given solution of HCl is taken in burette and run down drop wise in the conical flask, having 10ml MOH solution with few drops of phenolphthalein, till the end point is reached which is light pink. Three concordant readings are taken for volume of acid used. Observations

Obser va	<u>obset various</u>					
Sr. NO	Initial	Final reading	Volume of Acid			
	reading		used			
1	0	10	10 ml			
2	10	20	10 ml			
3	20	30	10 ml			
		Concordant I	Reading $(V_2) = 10 \text{ ml}$			
Calculations						
	Base		Acid			

Base		Acid
$\underline{M_1V_1}$	=	M_2V_2
n_1		n_2
<i>M</i> ₁ 10	=	0.1X10
1		1
M_1	=	0.1 X 10 X 1
		1 X 10
Molarity of MOH	=	0.1 M
Amount dissolve per	=	Molarity x Molecular.
Amount dissolve per litre of MOH	=	Molarity x Molecular. Wt.
-	=	-
litre of MOH		Wt.
litre of MOH	=	Wt. 0.1 x (M+16+1)
litre of MOH 4 <u>4</u>	=	Wt. 0.1 x (M+16+1)
litre of MOH 4 <u>4</u> <u>0.1</u>	=	Wt. 0.1 x (M+16+1) M+ 17
litre of MOH 4 <u>4</u> <u>0.1</u> 40	=	Wt. 0.1 x (M+16+1) M+ 17 M+ 17

Result The atomic weight of M is 23 g/mol. Which is considers as wt. of Na.

Practical Notes (Titration) F.Sc. 2nd year Chemistry

Page No. 7

Question: 06

The given solution contains 6g of a impure NaOH				
dissolve per litre. Determine the % age purity of				
the sample. You are provided with 0.1M HCl				
solution.				
Principle	It is an acid base titration.			
Standard Solution	0.1M HCl			
Chemical Equation	NaOH + HCl → NaCl +			
H ₂ O				
<u>Molar ratio</u>	$n_1 = 1, n_2 = 1$			
Indicator	Phenolphthalein			
End point	Light pink			
Procedure				

The given solution of HCl is taken in burette and run down drop wise in the conical flask, having 10ml impure NaOH solution with few drops of phenolphthalein, till the end point is reached which is light pink. Three concordant readings are taken for volume of acid used. <u>Observations</u>

	T *4* 1	T. 1 1.	X7 1 6 A 1 1
Sr. NO	Initial	Final reading	Volume of Acid
	reading		used
1	0	10	10 ml
2	10	20	10 ml
3	20	30	10 ml
		Concordant l	Reading $(V_2) = 10$ ml
<u>Calculat</u>	ions		
	Base		Acid
	M_1V_1	=	M_2V_2
	n_1		n_2
	<i>M</i> ₁ 10	=	0.1 <i>X</i> 10
	1		1
	M_1	=	0.1 X 10 X 1
			1 X 10
Mola	rity of NaOH	\sim	0.1 M
Amount	dissolve per	= Molar	rity x Molecular.
litre of N		Wt.	ity a molecului.
		=	0.1 x 40
		=	4g
6g of imp contains	oure sample	= 4	g of pure NaOH
1 a of im		_	4
contains	oure sample	=	$\frac{1}{6}$
100g of i	mpure sample	=	4
contains	mpare sumple	_	$\frac{4}{6}X$ 100
			66.6%

Result: % age purity of NaOH is 66.6%

UMAIR KHAN ACADEMY

www.youtube.com/umairkhanacademy www. UmairKhanAcademy.com

Question: 07	
<u>6.5 grams of an imp</u>	<u>ure sample of KMnO4 is</u>
<u>dissolves per litre. De</u>	<u>etermine the percentage of Mn</u>
present in the given s	<u>solution. 0.1M Mohar's satl</u>
<u>solution is given.</u>	
Principle	It is a Redox titration.
Standard Solution	0.1M Mohar's salt solution
Chemical Equation	
2KMnO ₄ + 10FeSO ₄ .(NI	$H_4)_2SO_4.6H_2O + 8H_2SO_4 \rightarrow$
$K_2SO_4 + 2MnSO_4 + 5Fe$	22 (SO4)3+10(NH4)2SO4 +68H2O
<u>Molar ratio</u>	
(KMnO4) n1 =2, (FeSO4	.(NH4)2SO4.6H2O) n2 =10
Indicator	KMnO4 Itself
End point	Light pink
Procedure	
The given sol	lution of KMnO4 is taken in

The given solution of KMnO4 is taken in burette and run down drop wise in the conical flask, having 10ml Mohar's salt solution and half test tube of dil. H₂SO4, till the end point is reached which is light pink. Three concordant readings are taken for volume of KMnO4 used. <u>Observations</u>

Sr. NO	Initial reading	Final reading	Volume of KMnO ₄ used
1	0	10	10 ml
2	10	20	10 ml
3	20	30	10 ml

```
Concordant Reading (V<sub>1</sub>) = 10 ml
```

Calculations		8. /
KMnO ₄		Mohar's Salt
M_1V_1	=	M_2V_2
n_1		n_2
<i>M</i> ₁ 10	=	0.1X10
2		10
\mathbf{M}_{1}	=	0.1 X 10 X 2
		10 X 10
Molarity of KMnO4	=	0.02 M
Amount of Mn dissolve per litre	=	Molarity x Atomic weight of Mn ⁺²
	=	0.02 x 55
	=	1.1g
% age of Mn in the given solution	=	$\frac{1.1}{6.5}X$ 100
	=	16.92%

Practical Notes (Titration) .8 F.Sc. 2nd year Chemistry

Page No. 8

Question: 08

The given solution contains 10 grams of mixture of ammonium oxalate and ammonium sulphate dissolved per litre. Find out the %age composition <u>of the mixture, you are provided with 0.02M</u> <u>KMnO4.</u> Principle It is a Redox titration. **Standard Solution** 0.02M KMnO4 solution **Chemical Equation** $2KMnO_4 + 8H_2SO_4 + 5(NH_4)_2C_2O_4.H_2O \rightarrow$ $K_2SO_4 + 2MnSO_4 + 5(NH_4)_2SO_4 + 10CO_2 + 13H_2O$ Molar ratio $(KMnO_4)$ n₂ =2, $(NH_4)_2C_2O_4H_2O$ n₁=05 Indicator **KMnO₄ Itself** End point Light pink Procedure

The given solution of KMnO₄ is taken in burette and run down drop wise in the conical flask, having 10ml of mixture solution and half test tube of dil. H₂SO₄+ heat to 60 ⁰C, till the end point is reached which is light pink. Three concordant readings are taken for volume of KMnO₄ used.

Observations

Obberva			
Sr. NO	Initial	Final reading	Volume of KMnO4
	reading		used
1	0	10	10 ml
2	10	20	10 ml
3	20	30	10 ml
		Concordant]	Reading (V ₂) = 10 ml
<u>Calculat</u>	ions		
(NH	[4)2C2O4.H2O	+ +	KMnO4
	M_1V_1	7	M_2V_2
	n_1		n_2
	<i>M</i> ₁ 10		0.02 X 10
	5		2
	M1	=	0.02 X 10 X 5
			2 X 10
Molarit	y of Amm.Oxala	te =	0.05 M
Amount dissolve	of Amm.Oxalate per litre	e = Mola	arity x Mol. weight
		=	0.05 x 142
		=	7.1g
10g of m Amm.Oxa	ixture contains date	=	7.1g
1g of mix Amm.Oxa	ature contains alate	=	7.1/10
100g of n Amm.Oxa	nixture contain date	s =	$\frac{7.1}{10}X$ 100
		=	71 g or 71%

UMAIR KHAN ACADEMY

www.youtube.com/umairkhanacademy www. UmairKhanAcademy.com

%age of Amm.Oxalate in the mixture	=	71%
%age of Amm sulphate in	_	100 51 200

%age of Amm.sulphate in = 100-71=29% the mixture

Question: 09

Determine the no. of molecules of water of crystallization in a given sample of oxalic acid by KMnO₄ titration. The amount of oxalic acid dissolved per litre is 6.3g. Principle It is a Redox titration. 0.02M KMnO4 solution **Standard Solution Chemical Equation** 2KMnO₄ + 3H₂SO₄ + 5(COOH)₂ .xH₂O → $K_{2}SO_{4} + 2MnSO_{4} + 10CO_{2} + (5x+8)H_{2}O_{3}$ Molar ratio $(KMnO_4) n_2 = 2, [(COOH)_2, xH_2O] n_1 = 05$ KMnO₄ Itself Indicator End point Light pink Procedure

The given solution of KMnO₄ is taken in burette and run down drop wise in the conical flask, having 10ml of oxalic acid solution and half test tube of dil. H₂SO₄+ heat to 60 ⁰C, till the end point is reached which is light pink. Three concordant readings are taken for volume of KMnO₄ used.

<u>Observations</u>					
Sr. NO	Initial reading	Final reading	Volume of KMnO4 used		
1	0	10	10 ml		
2	10	20	10 ml		
3	20	30	10 ml		

	Concor	dant Reading (V ₂) = 10 ml
Calculations		
(COOH) ₂		KMnO ₄
M_1V_1	=	M_2V_2
n_1		n_2
<i>M</i> ₁ 10	=	0.02 X 10
5		2
\mathbf{M}_{1}	=	0.02 X 10 X 5
		2 X 10
Molarity of Oxalic Acid	=	0.05 M
Amount of (COOH) ₂ dissolve per litre	=	Molarity x Mol. weight
6.3	=	0.05 x (90+18x)
126/0.05	=	90+18x
126	=	90+18x

=

=

18x

х

126-90

36/18

Practical Notes (Titration) F.Sc. 2nd year Chemistry

2

Page No. 9

x

=

Result: No. of molecules of crystallization in oxalic acid are two

Question: 10

The given solution contains 30g of partially oxidized FeSO₄.7H₂O dissolve per litre. Determine the %age oxidation of given Sample. It is a Redox titration. Principle 0.02M KMnO4 solution **Standard Solution Chemical Equation** 2KMnO₄ + 8H₂SO₄ + 10FeSO₄.7H₂O → $K_2SO_4 + 2MnSO_4 + 5Fe_2(SO_4)_3 + 78H_2O_4$ <u>Molar ratio</u> (KMnO₄) n₂=2, (FeSO₄.7H₂O) n₁=10 **KMnO4** Itself Indicator **End point** Light pink Procedure

The given solution of KMnO₄ is taken in burette and run down drop wise in the conical flask, having 10ml of partially oxidized ferrous sulphate solution and equal volume of dil. H₂SO₄, till the end point is reached which is light pink. Three concordant readings are taken for volume of KMnO4 used.

Observations	
Obsel valions	

<u>Observa</u>	<u>tions</u>			
Sr. NO	Initial	Final read	ling	Volume of KMnQ4
	reading			used
1	0	10		10 mł
2	10	20		10 ml
3	20	30		10 ml
		Concord	ant l	Reading (V ₂) = 10 ml
<u>Calculati</u>	ions			
Ferre	ous sulphate	\sim		KMnO ₄
	M_1V_1	⇒		M_2V_2
	$\overline{n_1}$			$\overline{n_2}$
	<i>M</i> ₁ 10	=		0.02 X 10
	10			2
	M ₁	=		0.02 X 10 X 10
				2 X 10
Molari	ty of Oxalic Acid	l =		0.1 M
Amount per litre	of FeSO4 dissol	ve =	Mola	arity x Mol. weight
		=		0.1 x 278
		=		27.8g
oxidized	unt of already l FeSO4.7H2O t Fe2(SO4)3	= 0		30 - 27.8= 2.2

UMAIR KHAN ACADEMY

www.youtube.com/umairkhanacademy www. UmairKhanAcademy.com

%age of Fe ₂ (SO ₄) ₃	=	$\frac{2.2}{30}X$ 100
	=	7.333%

Result: Percentage of oxidation is 7.333%

Ouestion: 11 <u>Determine the amount of iodine per litre in the</u> given sample solution. You are provided with 0.1M <u>Na₂S₂O₃ (hypo) solution.</u> Principle It is a Redox titration (Iodimetry) 0.1M Na₂S₂O₃ solution **Standard Solution Chemical Equation** $2Na_2S_2O_3 + I_2 \rightarrow Na_2S_2O_6 + 2NaI$ Molar ratio (Iodine) $n_1 = 1$, (Na₂S₂O₃) $n_2 = 2$ Indicator **Freshly prepared Starch Solution** just colourless End point Procedure

10ml of iodine solution is taken in conical flask and diluted it with a test tube full of water. Hypo (Na₂S₂O₃) solution is taken in burette and run drop wise in iodine solution until dark brownish colour of iodine changed to pale yellow colour and then about 2ml of freshly prepared starch solution as indicator is added which turned the solution colour into intensely blue. Now more hypo is run down drop wise until the solution changed to just colourless. Three concordant readings are taken for volume of hypo used

Observations

Sr. NO	Initial reading	Final reading	Volume of Hypo used
1	0	25	25 ml
2	10	35	25 ml
3	20	45	25 ml
		Concordant	Reading (V2) – 25 ml

Calculations

Concordant Reading $(V_2) = 25$ ml

Iodine		Hypo (Na ₂ S ₂ O ₃)
M_1V_1	=	M_2V_2
n_1		n_2
<i>M</i> ₁ 10	=	0.1 <i>X</i> 25
1		2
\mathbf{M}_{1}	=	0.1X25X1
		2 X 10
Molarity of iodine	=	0.125 M
Amount of Iodine dissolve per litre	=	Molarity x Mol. weight

=

0.125 x 254

Practical Notes (Titration) F.Sc. 2nd year Chemistry

Page No. 10

31.75g

Result: A/L of iodine is 31.75g

Ouestion: 12

Calculate the percentage purity of the given <u>impure sample of hypo, 80g of which are dissolved</u> per litre of the solution. You are provided with 0.05M iodine solution. Principle It is a Redox titration (Iodimetry) **Standard Solution** 0.05M iodine_solution **Chemical Equation** $2Na_2S_2O_3 + I_2 \rightarrow Na_2S_2O_6 + 2NaI$ Molar ratio $(Na_2S_2O_3) n_1 = 2$, (Iodine) $n_2 = 1$ Indicator **Freshly prepared Starch Solution** just colourless End point Procedure

=

10ml of iodine solution is taken in conical flask and diluted it with a test tube full of water. Hypo (Na₂S₂O₃) solution is taken in burette and run drop wise in iodine solution until dark brownish colour of iodine changed to pale yellow colour and then about 2ml of freshly prepared starch solution as indicator is added which turned the solution colour into intensely blue. Now more hypo is run down drop wise until the solution changed to just colourless. Three concordant readings are taken for volume of hypo used

Concordant Reading $(V_1) = 25$ ml Calculations Iodine Hypo (Na₂S₂O₃) M_1V_1 M_2V_2 n_2 n_1 0.05 X 10 *M*₁25 1 0.05 X 10 X 2 M_1 = 1 X 25 Molarity of hypo 0.25 M = Amount of hypo Molarity x Mol. = dissolve per litre weight 0.25 x 248 = 62g $\frac{62}{80}$ X 100 %age purity of hypo 77.5% =

Result: Percentage purity of Hypo is 77.5 %

UMAIR KHAN ACADEMY www.youtube.com/umairkhanacademy www. UmairKhanAcademy.com

Observations

Observa	10115		
Sr. NO	Initial	Final reading	Volume of Hypo used
	reading		
1	0	25	25 ml
2	10	35	25 ml
3	20	45	▲ 25 ml

Experiment 2

Experiment 1	
of Barium ions as Bariur	n chromate

Material Required

Ba2+ ions solution, K2CrO4 solution

Estimation

Apparatus

Beakern funnel tripod stand, burner wire gauze, oven, measuring cylinder, filter paper. etc.

Chemical Equation

 $Ba^{2+}K_2CrO_4 \longrightarrow BaCrO_4 + 2K^+$

Procedure

- i) Take 20 cm³ of Ba^{2+} ions solution in a beaker.
- ii) Now add K₂CrO₄ in excess in order to precipitate Ba²⁺ ions as BaCrO₄.
- iii) Boil this solution for 2-3 minutes to increase precipitate.
- iv) Filter the solution through a pre-weighed filter paper to get ppt.
- v) Wash the ppt. with hot water.
- vi) Dry the crystals and weigh them.

Observation & Calculation

Volume of Ba^{2^+} solution taken $= 20cm^3$ Mass of filter paper= 0.5gMass of filter paper + ppt.= 1.25gMass of BaCrO4 ppt.= 1.25 - 0.5 = 0.75gMolecular mass of BaCr4= 253.33 g/molAtomic mass of Ba²⁺= 137.33 g/mol

Thus

- 253.33g BaCrO4 contains Ba2+ ions
- 0.75 g BaCrO₄ contains Ba²⁺ ions

 20 cm^3 of the given solution contains Ba^{2+} ion 100 cm^3 of the given solution contains Ba^{2+} ions

Results

The given solution contains 2.03 g of Ba^{2+} ions in 100 cm³ of solution.

= 137.33g

0.406g

137.33×0.75

253.33

 $=\frac{0.406\times100}{20}=2.03$

= 0.406g

Detection of Elements in	an Oi	rganic	compound	(N,S,Hologen)
Prenaration of Lassaigne's s	olution	(LS)	•	

- i) Take a small piece of freshly cut sodium metal in a dry fusion tube.
- ii) Heat the fusion tube in the flame until sodium metal melts.
- iii) Add a small quantity of given organic compound in it.
- iv) Heat the mixture carefully till red hot.
- v) Put the fusion tube in a china dish containing about 10 cm³ of distilled water.
- vi) · Break the fusion tube and boil the solution for 5 minutes.
- vii) This solution is used for detection of elements (N,S & Halogens)
- viii) Filter the solutions. The filtrate called as Lassgne's solution or extract.

Test for nitrogen

Experiment	Observation	Inference
Sodium extract + few drops	i) Prussian blue colour or	Nitrogen is present.
of NaOH + freshly prepared	ppt. formed.	
$FeSO_4$ + boiled + cooled +	ii) Blood red colour	Both Nitrogen & sulphur
few drops of H ₂ SO ₄	appeared.	are present.
Test for Sulphur	•	i
Experiment	Observation	Inference
Sodium extract + Dil.	Black ppt. of lead sulphate	Sulphur is present.
CH ₃ COOH+ (CH ₃ COO) ₂ Pb	are formed.	
solution.	đ	N 1 2 2

Test for Halogens

Experiment	Observation	Inference
Sodium extract + few drops	i) White ppt. soluble in	Chlorine is present.
of Conc. HNO_3 + boiled +	excess of NH₄OH.	
$cooled + AgNO_3$ solution.	ii) Pale yellow ppt.	Bromine is present.
	sparingly soluble in NH ₄ OH	
	solution.	· · · · · · · · · · · · · · · · · · ·
	iii) Yellow ppt in insoluble	Iodine is present.
	in NH ₄ OH solution.	

Experiment 3

Identification of functional groups in a simple organic compounds

		ation of Carboxylic acid	
No.	Experiment	Observation	Inference
	Litmus Test Aq. Solution + blue litmus solution.	Blue litmus turned red.	May be an acidic compound (carboxylic acids).
2	<u>NaHCO3 Test</u> Aq. Solution + 5% NaHCO3 solution	A coloureless odourless gas evolved with effervescence, turned lime water milky (CO ₂ gas)	May be carboxylic acids.
3	Ester Formation Test Organic compound + Ethanol + few drops of Conc. H ₂ SO ₄	Fruity smell evolved due to ester formation	Carboxylic acid is confirmed.
4	FeCl ₃ Test Organic compound + few drops of neutral FeCl ₃	Red dish brown	Carboxylic acid is confirmed.
	i solution.		
Result	solution. The given organic co	mpound is a carboxylic acid.	
Result	The given organic co	ntification of Phenols	
	The given organic co ii) <u>Ide</u>	ompound is a carboxylic acid. ntification of Phenols Observation	Inference
<u>Result</u> <u>3r.No</u> 1	The given organic co ii) <u>Ide</u>	ntification of Phenols	Inference May be an acidic
3r.No	The given organic co ii) Ide Experiment Litmus Test Aq. Solution + blue	ntification of Phenols Observation	Inference May be an acidic compound (carboxylic
<u>3r.No</u> 1	The given organic colspan="2">ii) Ide Litmus Test Aq. Solution + blue litmus solution. NaHCO3 Test Organic compound +5%	ntification of Phenols Observation Blue litmus becomes red.	Inference May be an acidic compound (carboxylic acids, phenols). Carboxylic acids are absent.
<u>3r.No</u> 1 2	The given organic colspan="2">ii) Ide Experiment Litmus Test Aq. Solution + blue litmus solution. NaHCO3 Test Organic compound +5% NaHCO3 NaOH Test Organic compound +5%	ntification of Phenols Observation Blue litmus becomes red. No effervescences Organic compound is	Inference May be an acidic compound (carboxylic acids, phenols). Carboxylic acids are absent.
<u>3r.No</u> 1 2 3	The given organic colspan="2">ii) Ide Litmus Test Aq. Solution + blue litmus solution. NaHCO3 Test Organic compound +5% NaHCO3 NaOH Test Organic compound +5% NaOH. Solution. FeCl3 Test Organic compound + few drops of neutral FeCl3. Bromine Water Test Organic compound +	ntification of Phenols Observation Blue litmus becomes red. No effervescences Organic compound is soluble. Blue/purple colouration	Inference May be an acidic compound (carboxylic acids, phenols). Carboxylic acids are absent. Phenols are indicated.

Sr.No	Experiment	Observation	Inference
	Litmus Test Aq. Solution + blue litmus solution.	No change	Acidic compound absent.
1 .	NaHCO ₃ Test Organic compound +5% NaHCO ₃ solution	Insoluble	Carboxylic acids are absent.
2	NaOH Test Organic compound +5% NaOH. Solution.	No reaction occurred.	Phenols are absent.
3	2.4 Dinitrophenyl hydrazine test Organic compound + 2.4-DNPH	Yellow ppt. formed	Aldehyde or ketone is present.
4	Silver Mirror Test Organic compound + AgNO ₃ solution + 2-3 drops NH ₄ OH + warm on a water bath.	Silver mirror formed on the walls of the test tube.	Aldehyde is confirmed.
4	Fehling's solution test Organic compound + Fehling's solution + heat	Brick red ppt. of Cu_2O formed.	Aldehyde is confirmed.

PREPARATIONS OF

a) Aspirin

b) Iodoform

c) Glucosazone

d) Copper Ammine Complex

a) PREPARATION OF ASPIRIN Material required Material required Salicylic acid 5g, Acetyl chloride 4cm³, Pyridine 1cm³, Acetic acid 10cm³. Apparatus Apparatus Conical flask, measuring cylinder, beaker, water bath, funnel, filter paper etc. **Chemical equation** paper, wire gauze. COOH COOH **Chemical equation** CH1CH2OH + 412 + 3Na2CO3 -OH C-CH (Ethyl alcohol) (sodium carbonate) Pyridine + HCI Acetyl chloride Diagram Aspirin Salicylic acid Diagram water bath vire quaze Procedure i) bunsen burner Na₂CO₁ in it & shake well. Procedure, ii) Take a round bottom tlask. Add 5g of salicylic acid and h cm³ of pyridine i) it. in it. Then add 4 cm³ of acetyl chloride slowly with constant stirring. iii) 'ii)'' Maintain the temperature of mixture between \$0-60°C. iii) 149" "Now heat the flask in water bath for few minutes. iv) · Pour this mixture solution in about 200 cm³ of ice cold water & shake solution. V) well. v) White crystals of aspirin will be formed. vi) vi) Separate the crystals by filtration and wash them with cold water. vii) vii) Recrystallize crystals from equal volumes of water and acetic acid. viii) viii) Dry the crystals between folds of filter paper. ix)

Result


I prepared the pure orystals of aspirin and shown to the examiner:

b) PREPARATION OF IODOFORM

Ethyl alcohol z0cm³, sodium carbonate 20g, iodine 10g,

Round bottom flask, funnel tripod stand, measuring cylinder, water bath, filter

→ CHI₃ + 5Nal + 3CO₂ + 2H₂O + HCOONa (lodoform)



Take 100cm3 of distilled water in a round bottom flask. Add 20g of

Heat this solution on water bath upto 60°Cand add 20cm³ ethyl alcohol in

Add few crystals of solid iodine in above solution and shake well until the brown colour of iodine persist.

- At this stage yewllow crystals of iodoform begins to separate out from the
- Keep this solution for sometime .
- Separate these crystals by filtration.
- Recrystallize iodoform by dissolving in alcohol.
- Dry the crystals between folds of filter paper.

Result

I prepared the pure crystals of iodoform and shown to the examiner.

c) PREPARATION OF GLUCOSAZONE

Material required

Phenyl hydrazine 2g, glucose 1g, sodium acetate 3g, dilute acetic acid. Apparatus


Round bottom flask, water bath, funnel, filter paper, tripod stand, wire gauze etc. <u>Chemical equation</u>

H

$$C = O$$

 $HO - C - OH$
 $(CHOH)_3$
 $HO - C - OH$
 $(CHOH)_3$
 $H_2C - OH$
 $(Glucose)$
H
H
 $C = NNHC_6H_5$
 $C = NNHC_6H_5 + C_6H_5NH_2 + NH_3 + 2H_2O$
 $(CHOH)_3$
 $H_2C - OH$
 $H_2C - OH$

Diagram



Procedure

- i) Dissolve 1 g of glucose in 5cm³ water to form a solution.
- Take a round bottom flask and add 3g sodium acetate in 20cm³ of water.
 Add 2g of phenyl hydrazine into it to prepare a solution and shake well.
- -iii) Now add the glucose solution to the round bottom flask with constant shaking.
- iv) Heat the flask on a water bath. After about 20 to 30 minutes crystals of glucosazone start appearing.
- v) Cool the solution to get maximum crystals.
- vi) Filter the crystals and wash with dil. acetic acid.
- vii) Recrystallize the glucosazone with hot alcohol.
- viii) Dry the crystals between folds of filter paper.

Result

I prepared the pure crystals of glucosazone and shown to the examiner.

d) PREPARATION OF COPPER AMMINE COMPLEX

Material required Copper sulphate 5g, liquid ammonia 15cm³. Ethanol 40cm³. Apparatus Beaker, Measuring cylinder, stirrer, funnel, filter paper, watch glass c.c. Chemical equation CuSO₄ + 4NH₄OH \longrightarrow [Cu(NH₃)₄]SO₄ + 4H₂O Copper sulphate: [Cu(NH₃)₄]SO₄ + 4H₂O (Tetraammine copper (II) sulphat) (Copper ammine complex) Diagram

Copper mine complex-

Procedure

i) Take 20cm³ of water in a beaker and dissolve 5g CuSO₄ in it.

ii) Add few drops of Conc. H₂SO₄ in above solution to make it clear.

- iii) Then add liquid ammonia drop wise with constant stirring till an intense blue colour solution.
- iv) Now add 50 cm³ ot alcohol with constant stirring.
- v) Allow it to stand for sometime.
- vi) Long, thin blue needle shape crystals of copper ammine complex art formed.
- vii) Filter the crystal and wash with alcohol,
- viii) Dry the crystals between folds of filter paper.


Result

I prepared the pure crystals of copper ammine complex and shown to the examiner.

rv Tests CONCENTRATED GROUP Special group Palm Test Colour Flame test Filter ash Test SO; or PO, Dil Group Salt + conc H SO. Solid salt + Dil. H.SO. Blue Golden Pink Salt+H,O heat if necessary Cu" Or dil HCI Mg⁻¹ Yellow Na" Soluble For CH COOT reaction occur reaction occur dil. Group Is present above solution Vigorous Solid salt+oxalic NOTE NATURE OF GAS -, Slow Reaction reaction acid+drop of water BaCI soln. CH,COO, C,O, NO, Rub Effervescence CÍ Br Í Light Vinegar like smell Blue Odourless Violet White ppt . Note Nature of Gas Pink Al'or Note Nature of Gas Colourless CO, Gas Ma PO CO' or HCO' Above ppt. Rotten +Conc. HNO, Brown gas CH.COO Pungent egg smell smell,slow reaction Pugent NO. Yellow Red OR CH COOH Salt+H,O Violet Fumes Vinegar like smell Acetate Slow Reaction S-2 which turns K,Cr.O. smell give dense solution red Brown Fumes 1-0.S+dil, H.SO EXCESS Shake side of T.T Light USUMPLACE Palm test(+ve) C.O.2 1-0.S+ AgNO, +KMnO Soln. Green CH,COO' Green black Greenish Green white fumes (black ppt.) Slow reaction Colour Zn" Fe" Blue with NH,OH 1-0.S.+ FeCl,+ Boil reaction mixture, boil opt. Insoluble SO 2-0.S+CuSO, Light brown fumes discharged 2-0.S+FeSO,+ NI Cu* Salt Insoluble (Red colouration) Br' which increases (black ppt.) الك ما تم 1 . . 2-0.S+C,H,OH+H,SO, on addition of 1. 0. S. +AaNO, CO, present 3. 0. S+CaSO. SO; or S.O 1-0.S+AgNO, SE STERNE T.T. dil. H.SO. C1⁻¹ 1-0.S+AgNO, (No. Ppt.) (Fruity Smell) 8 yellow soln. (light yellow ppt.) paper pallet C,C . (yallow ppt.) 1-0.5+AgNO. 2. 0. S. +Lead acetate 2- 0.S+dil H,SO, confirmed blackening occur Salt soluble 1.S+CaCl, (white opt.) 2-0.5+dll H, SO,+ white ppt. (white ppt.) KMnQ+CCl+shake 2.S+dil. H,SQ,+KMnO, + MgSO, soln. KMnQ+CCI+shake For NH 2- chromyl 3. 0. S. +CaCl. (0. Layer Turns one drop + warm chloride test NO₃ 1- 0.S+Diphenyl amine (0. Layer turns Deep Orange) salt+Na_CO,+H_O (white ppt.) (salt solid + Dirty Check the violet) (blue colouration) Green Rub *Brick Red KCr.D.solid Colour of Greensolution of test 2- Ring Test Cr" Ppt. in + conc H.SO Ppt. on boiling KMnO. Smell of NH, Ca" (salt soln+EsSQ+H_SO, Conc along the side of tast tube brown Sn', Sn" Tube (orange red vapours) cold state CO decolourized HCO فرار من 1. 0.S+CdCl. -1. 0.S+CdCl,+boll ring is formed at the Soin, clear SO." ppt. Soluble (white ppt.) (white ppt.) tion of two solutio 1'0.S+BaCI. 2.0.S+BaCL . 2.0.S+BaCl,+boil PO (white ppt.) (white ppt) (white ppt.) · NH 2-0.S+ Lead actate Reddish 1 JAN SCALL (white ppt.) 1. 0. S. +AgNO, Pale Brown (yellow ppt.) Green الفاميد ادرتما توسلنيث كوبهت احتياطت جيك كري Fe 2. 0. S. +FeCl,(yellow ppL) Bat FILTER ASH TEST 3. 0. S. +Conc. HNO, salt + cobait nitrate solution +amm. Molybdate Solution + dil HCl (few drops) vellow ellow colour dip filter paper strip in S 0.* Salt. Soin+Co(NO,) above mixture. 1-0.S+ AgNO, (white , yellow , orange, black colour) Reddish Burn filter paper strip and 2-0.S+BaCl, (white ppt.) Pink Crimson 3-0.S+I Soln. (Decolorized) plet blue ppt: Co note the color of ash. Sr"

N-19

Page No. 15



Inference

		22	4 -		4 1
	1		Bead or Scales	Incrustation	
		1	White grey bead, marks paper.	Yellow incrustation.	Pb
		2	White malleable bead, does not	No incrustation	Ag
			mark paper.		
	Q.7	Exp	lain chemistry of Charcoal Cavit	v Test?	
	Ans:	Na ₂	CO3 converts a metallic salt into a	carbonate which decomp	oses to its oxide
acid with a	on heating. The oxide in certain cases is further reduced with the carbon of				
i acid with a		char	rcoal block to metal which appears i	n the form of bead or scal	.es.
			$Pb(NO_3)_2 + Na_2CO_3 \longrightarrow 21$		
			PbCO ₃ → PbO + CO ₂	(Yellow incurstation	ı) ·
			(Bead)		
ical. It is			$PbO + C \longrightarrow CO + Pb$		
n-metals e.g.		In c	certain cases metal oxides are no	t reduced to the metal	or only partially

Observations

$$CdCl_2 + Na_2CO_3 \longrightarrow CdCO_3 + 2NaCl$$

 $CdCO_3 \longrightarrow CO_2 + CdO$ (Brown incrustation).

What is the difference between a bead and an incrustation?

Ans: Bead is a spherical shaped metallic mass. Incrustation is a layer of metallic oxides formed around the cavity.

 $Pb(NO_3)_2 + Na_2CO_3 \longrightarrow 2NaNO_3 + PbCO_3$

$$PbCO_3 \longrightarrow PbO + CO_2$$
 (Yellow incurstation)
(White grey head)

Q.9 What is Filter Ash Test

Q.8

Ans: Crystals of salt are dissolved in equal amount of cobalt nitrate solution. A filter paper strip is dipped in this solution. Dry it, burn on flame and colour of ash is noted.

Observation	Result	
1) Pink ash	Mg Salt	
2) Green ash	Zn salts	
3) Blue ash	$Al^{3+}+ Or PO_{A}^{-3}$	
4) Dirty green ash or Bluish green	Sn salts	

Q.10 How would you perform flame Test?

Ans: Make thick paste of salt with drops of Conc. HCl in China dish (salt is converted to chlorides which are more volatile compound). With the help of platinium wire or scapula, this paste is placed in oxidizing flame

Selective Short questions

Q.1 What is salt?

Ans: <u>Salt</u>

Salt is a compound which is formed by complete neutralization of an acid with a base. NaOH + HCl \longrightarrow NaCl + H₂O

e.g. NaCl, NaHCO3, KNO3 etc.

Q.2 What is acidic and basic radical?

Ans: Acidic radicals

That part of the salt which is derived from the acid is called acid radical. It is negatively charged portion of the salt. It is usually composed of a non-metals e.g. $CO_3^{2-}, SO_4^{2-}, B\overline{r}, C\overline{l}$ etc.

Basic Radical

That part of the salt which is derived from the base. It has positive charge. It is usually a metallic ion e.g. Na^+ , K^+ , Ca^{2+} , Al^{3+} etc. or NH_4^+ .

Q.3 What is the difference between the detection of basic radicals and acid radicals?

Ans: The basic radicals are detected and confirmed by the formation of colour or precipitate, while the acid radicals are generally detected by the confirmation of their volatile product.

Q.4 What do you mean by an "ion".

Ans: An ion is an atom or a group of atoms having positive or negative charge. e.g. CO_{1}^{2-} , NO_{1}^{1-} , CI^{-} , SO_{4}^{2-} , Na^{+} , K^{+} etc.

Q.5 What is chemical test?

Ans: A chemical reaction performed on the salt or on its solution, to check presence or absence of some radical is called chemical test.

Q6 How would you perform charcoal cavity test?

Ans: Dip up a shallow cavity in a charcoal block. Mix a little salt with anhydrous Na₂CO₃ and fill the charcoal cavity with it. Put a drop of water on the salt to make a paste. Heat it strongly in reducing flame of Bunsen burner with the help of blow pipe. Examine the bead and colour of the incrustation or scales.

N-21

Page	No.	18
------	-----	----

Observation	Result	
1) Bluish green	Copper Cu ²⁺	
2) Yellowish green or grass green	Barium Ba ²⁺	
3) Crimson	Strontium Sr2+	
4) Brick red	Calcium Ca ²⁺	
5) Violet	Potassium K ⁺	
6) Golden yellow	Sodium Na ⁺ ·	

Q.11 What is the purpose of using blue glass in the flame test?

- Ans: Blue glass absorbs a part of light. Therefore, blue glass produces a characteristic change in the colour of the flame of basic radicals having similar flames. For example Ca and Sr have brick red and crimson colour and they appear green and purple through blue glass.
- Q.12 Explain Chemistry of Borax Bead Test?

Ans:
$$Na_2B_4O_7 \xrightarrow{} 2NaBO_2 + B_2O_3$$

(Volatile) $CuSO_4 \longrightarrow CuO + SO_3$ $B_{2}O_{3} + CuO \longrightarrow Cu(BO_{2})_{2}$ (blue bead) $2Cu(BO_{2})_{2} + C \longrightarrow 2CuBO_{2} + B_{2}O_{3} + CO$ (Colourless) $2CuBO_2 + C \longrightarrow 2Cu + B_2O_3 + CO$

(Reddish brown)

Q.13 Can you perform borax bead test for colourless salts?

- No, this test is only performed for coloured salts because none of the colourless Ans: salt gives a characteristic bead with borax beat test.
- 0.14 Why is platinum wire used in the borax bead test?
- Only platinum wire is used in the borax bead test because it has high melting Ans: point and imparts no colour to the flame.
- Q.15 Can we make a paste of the salt in Conc. H2SO, instead of Conc. HCl?
- Ans: No. if a paste of the given salt is prepared with Conc. H2SO4, it will change to sulphate which are not as volatile as chlorides. Therefore characteristic colour will not be imparted to the flame easily.
- Q.16 What do you mean by a wet test?
- Ans: It is attest which is applied to the salt after forming its solution. e.g. confirmatory tests of acid & basic radicals.
- Q.17 Explain preparation of Neutral Solution (N.S.)?
- Ans: Few gram of solid salt + Double quantity of solid Na₂CO₃ + 5cm³ H₂O + Boil for 10-15 minutes then centrifugate it, discard the residue. Centrifugate is basic due

to addition of Na2CO3 Then cool the centrifugate + CH3COOH with stirring until evolution of CO2 gas ceases (slightly acidic) + NH4OH with stirring. Boil the soln. to expel the excess of NH₃ gas, then neutral solution is obtained which is used for confirmatory test.

- 0.18 Why do we heat the salt with strong solution of Na₂CO₃ for the preparation of solution of acid radicals?
- Ans: We heat the salt with strong solution of Na₂CO₃ so that double decomposition reaction may take place. The anion remain in the solution while cations of heavy metals are replaced by sodium and are precipitated as neutral or basic carbonates and in some cases hydroxides. For example,
 - $BaSO_1 + Na_2CO_1 \longrightarrow$ $BaCO_1 \downarrow + Na_2SO_1$ i)
 - $2CuSO_4 + 2Na_2CO_3 + H_2O \longrightarrow (CuOH)_2CO_3\downarrow + 2Na_2SO_4 + CO_2\uparrow$ ii)
 - $2AICl_3 + 3Na_2CO_3 + 3H_2O \longrightarrow 2AI(OH)_3 \downarrow + 6NaCl + 3CO_2\uparrow$ iii)
- 2.19 Name the acid radicals of dilute acid group.
- Ans: (is) Carbonate (CO_1^{2-}) (iii) Sulphide (S²⁻)
 - (iv) Sulphite $(SO_3^{2^-})$

(ii) Bicarbonate (HCO_3^{I-})

- (v) Thiosulphate $(S_1O_1^{2-})$ (vi) Nitrite (NO; -)
- Q.20 Name some insoluble carbonates.
- All carbonates are insoluble in water except those of sodium, potassium and Ans: ammonium. For example, CaCO₃, SrCO₃, BaCO₃ etc.
- 0.21 Why lime water turns milky on passing CO2 through it?
- Ans: When CO2 is passed through lime water Ca(OH)2, insoluble CaCO3 is produced as white milky suspensions.
 - $Ca(OH)_2 + CO_2 \longrightarrow CaCO_3 + H_2O_3$
- Q.22 If an excess of CO₂ gas is passed through the above milky precipitate, the milky precipitate disappears and a clear solution is formed why?
- Ans: On passing excess of CO2, the precipitate gets dissolved due to the formation of Ca(HCO₁)₂
 - $CaCO_3 + H_2O + CO_2 \longrightarrow Ca(HCO_1)_2$
- Is there any gas like CO2, which turns lime water milky? 0.23
- Ans: Yes, SO2, turns lime water milky due to the formation of insoluble CaSO1. $Ca(OH)_2 + SO_2 \longrightarrow CaSO_3 + H_2O$
- 0.24 How will you distinguish between a carbonate and a bicarbonate?
- When MgSO4 is added to O.S. formation of a white ppt. in the cold state indicates Ans: carbonate. But the formation of a white ppt. on heating indicates the presence of bicarbonate.

 $CO_3^{2-} + MgSO_4 \xrightarrow{\ln \infty M} MgCO_3$ (White ppt.) $HCO_3^{l-} + MgSO_4 \longrightarrow Mg(HCO_2)_2 \longrightarrow MgCO_3 + CO_2 + H_2O$ (White ppt.)

- Why a bicarbonate forms white ppt. on heating with MgSO4? Q.25
- Because all bicarbonates are soluble in water and on heating they are changed to Ans: insoluble carbonates.

 $MgSO_4 + 2NaHCO_3 \rightarrow Na_2SO_4 + Mg(HCO_3)_2$ $Mg(HCO_3) \longrightarrow MgCO_3 + CO_2 + H_2O$ White ppt.

How will differentiate between SO_1^{2-} and $S_1O_2^{2-}$? 0.26

Ans: Add dil. H₂SO₄ to the salt. A colourless gas with burning sulphur smell which turns K2Cr2O7 solution green indicates both sulphite and thiosulphate. Note the colour of the contents of the tube.

Yellow contents = Thiosulphate

Transparent = Sulphite

- 0.27 What is palm test?
- When solid salt containing acetate radical is mixed with solid oxalic acid in the Ans: presence of a drop of water, smell of acetic acid or vinegar comes out. This test usually performed on palm of the hand. So it is called palm test.
- What is the fruity smell test? 0.28
- When the solution of an acetate salt is heated in the presence of ethyl alcohol and Ans: Conc. sulphuric acid, fruity smell comes out.

 $CH_3COO^{1-} + C_2H_5OH + Conc, H_2SO_4 \longrightarrow CH_3COOC_2H_5 + H_2O_{1-}$

(Fruity smell)

- Does nitrite gives "ring test"? 0.29
- Nitrite also gives ring test similar to that of nitrates but the acid used should be Ans: dil. acetic acid or dil. H2SO4 instead of Conc. H2SO4.

3HNO, H₂O + HNO₃ + NO

FeSO4 + NO _____ [Fe(NO)]SO4

Nitroso-ferrous sulphate (ring)

- Q.30 Can we use Conc. HCl in place of Conc. H_SO4 in concentrated acid group.
- Ans: No, we cannot use Conc. HCI in Conc. acid group because this group contains CI ion so it will create confusion because we are introducing Cl ions ourself in the form of Conc. HCl.
- Q.31 What is chromyl chloride test?
- It is confirmatory test for Cl and is performed as. Take small quantity of solid Ans: salt and solid K₂Cr₂O₇ in a test tube. Add Conc. H₂SO₄ in the test tube and heat the mixture. Pass the gas through NaOH solution and then add acctic acid and lead acetate. Yellow ppt. of PbCrO4 is obtained.

 $K_2Cr_2O_7 + 4NaCl + 6H_2SO_4 - \frac{\Lambda}{Heat}$ KHSO₄ + 4NaHSO₄ + 5H₂O + 2CrO₂Cl₂ \uparrow (Chromyl chloride) CrO₂Cl₂ + 4NaOH Na2CrO4 + 2NaCl + 2H2O $Na_2CrO_4 + (CH_3OO)_2Pb$ PbCrO₄ + 2CH₃COONa (Yellow ppt.) How the reddish brown vapours of Br2 and NO2 are distinguished? 0.32

Bromides give the brown vapour of bromine while nitrate gives brown vapours Ans: of NO2. If the addition of Cu-turnings, increases the intensity of brown fumes, then nitrate is indicated and if there is no change in vapour's intensity on the addition of Cu-turnings, then bromide is present.

- By what methods SO²⁻ and PO³⁻ are confirmed? 0.33
- The confirmation of these two radicals differs from others in the fact that they are Ans: confirmed by the method of precipitation and not by volatile products.
- What happens when silver nitrate solution is added to a thiosulphate? 0.34
- Ans: A white precipitate of silver thiosulphate is formed which immediately changes its colour to yellow, organs, brown and finally black.

 $S_2O_1^{2-} + \text{AgNO}_3 \longrightarrow \text{Ag}_2S_2O_3$

(White coloured)

 $Ag_2S_2O_3 + S_2O_3^{2-} \longrightarrow [Ag(S_2O_3)]^{-3} \longrightarrow Ag_2S_2O_3 + S_2O_3^{2-} \longrightarrow Ag_2O_3^{2-} \longrightarrow Ag_2O_3^{2-} \longrightarrow Ag_2O_3^{2-} \longrightarrow Ag_2O_3^{2-} \longrightarrow Ag_2O_3^{2-} \longrightarrow Ag_2O$

(Black coloured ppt.)

- Q.35 Name the acid radicals of Conc. Sulphuric acid group.
- Ans: $Cl^{1-}, Br^{1-}, I^{1-}, NO_1^{1-}, CH_1, COO^{1-}, C_2O_2^{1-}$
- Q.36 How will you distinguish between a thiosulphate and sulphite by the help of dil. H2SO4 only?
- Dilute sulphuric acid will produce SO₂ gas only and the solution in the test tube Ans: will remain transparent . But in case of thiosulphate the contents of the test tube will become yellow due to formation of a ppt. of sulphur along with evolution of SO₂ gas.

$$SO_3^{2^-} + H_2SO_4 (dil.) \longrightarrow SO_4^{2^-} + SO_2\uparrow + H_2O$$

$$S_2O_3^- + H_2SO_4(dil.) \longrightarrow SO_4^- + SO_2 \uparrow + S\downarrow + H_2O$$

0.37 What is common ion effect?

Ans: Common ion effect is a phenomenon in which the degree of dissociation of a weak electrolyte is suppressed by the addition of another electrolyte having common ion.

> $H_2S \implies 2H + S^{2-}$ $HCI \longrightarrow H^+ + CI^-$

- What is Original solution:(O.S.) 15
- Ans: A transparent, homogeneous solution of salt in a suitable solvent. This is used for analysis of avid and basic radicals by wet tests. If dissolved in water, also called aqueous solution (Aq.S.)
- 0.39 Sometimes we get a white precipitate on diluting a salt solution prepared in Conc. HCl. Explain why?
- Ans: If the salt is soluble in Conc. HCl but forms white ppt. on dilution, it may be due to formation of oxychlorides of Bi¹⁺ or Sb¹⁺ e.g.

BiCl1 + H2O ----- BiOCI + 2HCI

White ppt.

The white ppt, can be redissolved by adding few drops of Conc. IICI.

BiOC1 + 2HC1 ----- BiCl1 + H2O

- 0.40 Why is AgNO₃ solution kept in coloured bottles?
- Ans: Since AgNO1 decomposes into its oxides in direct sunlight, therefore, it is kept in coloured bottles to prevent its decomposition.
- 0.41 How will you perform the lake test for Al?
- Ans: Jo O.S. solution, add few drops of litmus solution and acidify it with dil. HCl. Then add NH4OH. Allow it to stand for a while. A blue ppt. floats in a colourless solution.
- How will you differentiate between Fe2' and Fe3'? 0.42

Ans:

Fe ²⁺	(i) (ii)	O.S. + NaOH soln. O.S. + potassium ferricyanide K ₃ [Fe(CN) ₆].	Yellowish green ppt. Dark blue ppt.	Fe ^{2*} is confirmed Fe ^{2*} is confirmed
Fe ³⁺	(i) (ii)	O.S. + NaOH soln. O.S.+Pot.ferricyanide + K ₃ [Fe(CN) ₆].	Blue ppt. Deep blue ppt.	Fe ³⁺ is confirmed Fe ³⁺ is confirmed

- Q.43 Are the basic radicals metals only?
- Ans: The majority of the basic radicals are metals but non-metallic radicals like Ammonium, are also basic radicals.
- Q.44 Why some radicals are coloured and some are not?
- Ans: Metal of the d-block series, which have unpaired electrons in their valence shell, show colour due to transition of these electron in the d-orbitals. Anions of transition metal complexes like permanganate ion or the chromate ion show colour due to the movement of odd electron in the molecule.
- Q.45 What are deliquescent substances?
- Ans: These are the substances which absorb moisture from the air. For example, NaOH, P2O5 etc.

- Name the basic radicals of Group 1" ? How the basic radicals of group 1st are 0.46 precipitated?
- Ans: Basic radicals of group I are Pb2', Hg1' Ag1'. The basic radicals of group 1st are precipitated by dilute HCl.
 - Basic radicals of group I are precipitated as their chlorides, e.g. AgCl, HgCl, PbCl2.
- Q.47 Can we use Dil. H2SO4 in place of Dil. HCl as group reagent in first group of basic radicals?
- Ans: No, because in that case basic radicals of other groups will also be precipitated as sulphates in the first group, e.g. BaSO4, SrSO4 etc.
- 0.48 Name basic radicals of group 2ªd. What is the group reagent for 2ªd group?
- Ans: 2" A Pb4", Hg2", Bi3", Cd2", Cu2". 2" B Sn2+, Sn4+, As3+, Sb3+. H₃S in the presence of Dil. HCl.
- 0.49 Why yellow ammonium Sulphide is used?
- Ans: It is used in the 2rd group of basic radicals to distinguish between basic radicals of HA and IIB groups.
- 0.50. Why do we add dil. HCl before passing H2S gas in the 2nd group?
- Ans: H₂S is a weak electrolyte and cannot ionize completely in the presence of common H' ions from HCl which is a strong electrolyte.
 - $H_{rS} \longrightarrow 2H + S^{2-}$

HCI ==== H' + CI

This low concentration of sulphide ions would precipitate basic radicals of 2nd group only, whereas it would be insufficient to precipitate radicals of other groups because their sulphides are comparatively soluble.

- Q.51 .Sometime a turbid solution is obtained when the mixture is dissolved in water, what is that due to?
- Ans: Due to Bi and Sb salts which hydrolyse to form BiOCI and SbOCI. It may be avoided by adding few drops of dil. HNO3, which will prevent hydrolysis.
- 0.52 Name basic radicals of 3rd group. What is group reagent of 3rd group.
- Ans: Fe2+, Fe3+, Al3+, Cr3+. NH4Cl(1) + boil + cool + NH4OH.
- 0.53 Why NH4Cl is added along with NH4OH?
- The cations of 3rd group are precipitated as insoluble hydroxides, if NH4Cl is not Ans: added other cations may also precipitate as hydroxides. In the presence of NH4Cl, the ionization of NILOH is suppressed so that the -OH ions in the solution are just sufficient to precipitate the radicals of the 3rd group only and not others.

$$H_4Cl \implies NH_4^* + Cl^-$$

NH4OH NH: + OH-

Q.64 Name the basic radicals of 6th group. What is the group reagent of 6th group? 6.54, What is green vitriole? Ans: Mg^{2+} , Na^+ , K^+ and NH_4^+ FeSO₄.7H₂O is called green vitriole. Ans: There is no common group reagent for the basic radicals of 6^{th} group. Hence they What is alum? 0.55 are confirmed individually except Mg2+. Mg2+ is precipitated by the above Ans: It is a double salt of potassium an aluminium having the formula K2SO4 Al2(SO4)3.24H2O. solution of 5th group + (NH4)3PO4 Q.56 Why is it necessary to boil off H₂S from the filtrate of 2nd group before the Q.65 Name some colour salts of Na and K. addition of Conc. HNO1? $K_2CrO_4 = Yellow$ Ans: $Na_2CrO_4 = Yellow$ Ans: It is essential to boil off H2S gas, otherwise it will be oxidized, at least partly, to $K_2Cr_2O_7 = Orange red.$ $KMnO_4 = Pink$ H2SO4. 0.66 What is the formula of diphenyl amine & dimethylglyoxime. $H_2S + 8HNO_3 \longrightarrow H_2SO_4 + 4H_2O + 8NO_2$ Ans: (C₆H₅)₂NH or If the mixture contains Ba and Sr, they will form insoluble sulphates and get precipitated in the 3rd group. $CH_3 - C = NOH$ (ii) Q.57 Name the basic radicals of 4th group. What is the group reagent for 4th group? (i) $CH_3 - \dot{C} = NOH$ Ans: 4th group of basic radical contains Co²⁺, Ni²⁺, Zn²⁺ and Mn²⁺. The group reagent is $NH_4Cl + boil + cool + NH_4OH + H_2S$. 0.58 How the basic radicals of 4th group are precipitated? What is the formula of (i) Nessler's reagent (ii) Picric acid, (iii) sodium 0.67 Ans: They are precipitated as their sulphides. cobaltinitrite, (iv) Potassium pyroantimonate? For example MnS = Pink or flesh colour ppt. Ans: i) Nessler's reagent: K₂Hgl₄ NiS = Black ppt. ZnS = Grayish white ppt. ii) Picric acid. OH CoS = Black ppt. 0.59 What are the colour of Cu, Fe, Co, Ni, Cr and Mn beads? Fe = reddish brown bead i) Cu = blue bead(ii) Ans: Ni = Brown bead Co = blue bead (iv) iii) Cr = Deep green bead. Mn = violet bead (vi) NO₂ (v) iiii) Sodium cabltintrite Na₃[Co(NO₂)₆] 0.60 What are the basic radicals of group 5th? What is the group reagent of 5 iv) K₂H₂Sb₂O₇ is potassium pyroantimonate. Ba^{2+} , Ca^{2+} and Sr^{2+} . Ans: (NH₄)₂ CO₃ in the presence of NH₄Cl and NH₄OH. Q.61 How the basic radicals of 5th group are precipitated? Ans: Basic radicals of 5th group are precipitated as their carbonates. CaCO₃, SrCO₃, BaCO₃. 0.62 Why is it necessary to add NH4OH before the addition of (NH4)2CO3 for the precipitation of group 5th radicals? (NH4)2CO3 reagent used in the laboratories usually contains a high percentage of Ans: NH4HCO3 which will give a little precipitate with the metals of group 5th if NH4OH is not added. \rightarrow (NH₄)₂CO₃ + H₂O NH4HCO3 + NH4OH Q.63 Sometimes it is expected that calcium does not get precipitated in group 5th. What is the reason? Ans: This happens because Ca forms Ca(HCO₃)₂ instead of CaCO₃ on the addition of (NH₄)₂CO₃ solution. Ca(HCO₃)₂ is soluble and passes into the filtrate of V group.